



**INTERNATIONAL JOURNAL OF
PHARMACEUTICAL SCIENCES**
[ISSN: 0975-4725; CODEN(USA): IJPS00]
Journal Homepage: <https://www.ijpsjournal.com>



Review Article

Instrumental Approaches in Neuropharmacological Screening: A Review of Insights from Electroconvulsimeter and Actophotometer Studies

M. Nithya*, G. M. Sivakumar, G. Kiruthika, A. Kokila, A. Mohanasundaram, S. Prasanna, N. Rithikadevi, V. Selvaraja, S. R. Shibi, T. Thirisha

JKK Munirajah Institute of Health Sciences College of Pharmacy, T. N. Palayam, Erode, Tamil Nadu, India 638506

ARTICLE INFO

Published: 25 Aug 2025

Keywords:

Electroconvulsimeter, Actophotometer, CNS drugs, Tonic-clonic seizures and Rodents

DOI:

10.5281/zenodo.16940508

ABSTRACT

The Electroconvulsimeter and Actophotometer are key experimental instruments in neuropharmacological research, widely employed for screening and evaluating the effects of central nervous system (CNS)-active agents. The Electroconvulsimeter functions by inducing controlled maximal electroshock seizures (MES) in rodents, providing a reliable and reproducible model for assessing the anticonvulsant efficacy of drugs, particularly those targeting generalized tonic-clonic seizures. The actophotometer assesses locomotor activity through infrared beam breaks produced by animal movement. This technique is valuable for evaluating the stimulant or depressant actions of CNS drugs, offering a non-invasive and quantitative method of behavioral analysis. Together, these instruments serve as indispensable tools in preclinical pharmacology, contributing to the discovery, characterization, and optimization of therapeutic compounds by facilitating objective assessment of their anticonvulsant, sedative, stimulant, and neuroprotective properties.

INTRODUCTION

Electroconvulsimeter¹⁻¹⁰:

The Electroconvulsimeter is a specialized scientific instrument used in neuropharmacological research to study seizure

activity and evaluate the anticonvulsant potential of pharmaceutical compounds. This instrument delivers controlled electrical stimuli to laboratory animals, typically mice or rats, to induce convulsions, thereby mimicking epileptic conditions in humans. The induced seizures allow researchers to assess the protective efficacy of

*Corresponding Author: M. Nithya

Address: JKK Munirajah Institute of Health Sciences College of Pharmacy, T. N. Palayam, Erode, Tamil Nadu, India 638506

Email ✉: nithi.m96@gmail.com

Relevant conflicts of interest/financial disclosures: The authors declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.



antiepileptic drugs by observing the alteration or suppression of convulsive phases such as tonic and clonic seizures. The Maximal Electroshock Seizure (MES) model, facilitated through the Electroconvulsimeter, is one of the most widely accepted and reproducible models for screening antiepileptic agents, especially those effective against generalized tonic-clonic seizures.

The Electroconvulsimeter remains an indispensable tool in early-phase drug discovery and in understanding the neurobiological basis of epilepsy. Its relevance is further underscored by the high translational value of the MES model, as several clinically used antiepileptic drugs such as phenytoin and carbamazepine were originally identified using this method. The Electroconvulsimeter continues to play a central role in pharmacological screening due to its reliability, simplicity, and cost-effectiveness. The Electroconvulsimeter operates on the principle of electrically induced seizures in experimental animals to evaluate the anticonvulsant potential of drugs. When a precisely controlled alternating current (AC) of defined intensity, frequency, and duration is passed through the brain via corneal or ear-clip electrodes, it causes widespread depolarization of neuronal membranes and synchronous firing of neurons. This results in a generalized tonic-clonic seizure pattern known as the Maximal Electroshock Seizure (MES). The MES consists of several phases, the most characteristic being Tonic Hind-Limb Extension (THLE), where the hind limbs extend at an angle of $\geq 180^\circ$ to the body axis. Anticonvulsant drugs that are effective against generalized tonic-clonic epilepsy are able to prevent or shorten the THLE phase, and this forms the basis for their screening using the Electroconvulsimeter.

Major Components of Electroconvulsimeter:

The major components of an Electroconvulsimeter are as follows:

1. **Stimulus Generator** is the primary unit of the Electroconvulsimeter. It is responsible for producing an electrical stimulus with controlled parameters such as: Current intensity (commonly 50–150 mA), Frequency (typically 50–60 Hz), and duration (0.2 to 0.5 seconds). These parameters can be finely tuned to ensure consistency across experiments. The generator typically uses square wave pulses, as they are more effective in eliciting maximal electroshock seizures (MES).
2. **Electrodes** are used to deliver the electrical stimulus to the animal. They are generally made from stainless steel or other biocompatible metals and are of two main types: Ear clip electrodes – attached to the pinnae of the animal. Corneal electrodes applied to the cornea after moistening with saline to reduce resistance. The appropriate electrode is selected based on the species and responsiveness of the subject.
3. **Current/Voltage Regulator** is primarily used as component ensures that the voltage and current delivered are within safe and precise limits. It allows researchers to adjust the output depending on the experimental requirements. Safety cut-offs are often built-in to prevent damage to animal tissues.
4. **Timer and Control Switches** may be a digital or analogue timer and is crucial for maintaining uniform stimulation periods across tests. Control switches allow for starting, stopping, or resetting the device.
5. **Display and Indicators** provides real-time feedback on key parameters such as current, voltage, frequency, and stimulus duration. This ensures accurate monitoring during the procedure.



6. Animal Restraint Device is not part of the electrical circuitry, animal restraint systems (e.g., small cages or holders) are essential for: Preventing excessive movement, ensuring accurate electrode placement, Enhancing animal safety and experimental reproducibility.

Standard Operating Procedures:

1. Divide the experimental animal into 3 groups (vehicle control standard drug).
2. Weight each experimental animal and mark them accordingly.
3. Administer the drug (standard drug and test drug) and normal saline vehicle control 60 min prior if given oral and 30 min prior if given I.P.
4. Keep the animal in the restrainer (for the ease for electrode fixation) or hold the animal properly.
5. Ensure that the place of contact for electrodes (ear electrode) is clean and apply electrode gel for better electrical conduction.
6. Place the electrode (ear electrode) on both ear pinna.
7. Set the desired electroshock parameters in Electroconvulsimeter.
8. Switch on the main power and press the start button.
9. Turn on the stop watch as soon as the start button is pressed and note the time spend in different phases of MES convulsion and note it down.
10. The Electroconvulsimeter will auto-off once the set time is completed.
11. Place the animal back into cage once it has recovered from the convulsion.
12. Repeat the same steps for the other groups.

Animals Used:

1. Rats (Albino Sprague–Dawley or Wistar)

Commonly used for Maximal Electroshock Seizure (MES) testing. Adult rats weighing 150–250 g are preferred, Stimulation parameters: 150 mA, 50–60 Hz, 0.2s duration. Electrodes: corneal or ear-clip electrodes. It is primarily used to induce tonic hind-limb extension for evaluating anticonvulsant drugs.

2. Mice (Albino strain)

Typical age: 6–8 weeks, weight 20–30 g, MES protocol: 50 mA, 60 Hz, 0.2 s duration. Electrodes: corneal type, with saline or tetracaine applied to improve contact. Purpose: Screening of potential anticonvulsant agents by observing seizure prevention or suppression.

Applications:

1. Screening of anticonvulsant drugs using the Maximal Electroshock Seizure (MES) model.
2. Evaluation of CNS depressant or stimulant effects of various pharmacological agents.
3. Determination of seizure threshold and convulsion severity in experimental animals.
4. Differentiation between types of anticonvulsant drug mechanisms (e.g., sodium channel blockers vs. GABA enhancers).
5. Comparative analysis of drug efficacy in suppressing electrically-induced seizures.
6. Assessment of the neuroprotective potential of new synthetic or herbal compounds.
7. Pharmacological interaction studies involving anticonvulsant agents.
8. Standardization and validation of seizure models for research purposes.
9. Educational demonstrations of seizure induction and drug action in pharmacology labs.
10. Evaluation of toxicity or side effects of CNS-active compounds through seizure susceptibility.





Fig.1 Diagramatic representation of Electroconvulsometer

Actophotometer¹¹⁻¹⁶

Actophotometer is an instrument used to measure locomotor (activity) behaviour in animals. It is commonly used in CNS drug screening to assess stimulant or depressant effects. The instrument detects animal movement using infrared light beams and photocells. Each time the animal crosses a beam, it interrupts the light, generating a count. The total count over a fixed period reflects the level of motor activity. Increased activity suggests CNS stimulation (e.g., with caffeine), while reduced activity suggests CNS depression (e.g., with diazepam). It is useful in testing neuroactive drugs, sedatives, anxiolytics, and stimulants.

The instrument can be digital or manual, with circular or square enclosures. Both rats and mice can be tested using different models. It provides a non-invasive and reliable method to quantify spontaneous movement. It helps to quantify drug induced changes in animal motor behaviour objectively. It is widely used in preclinical pharmacology and behavioral neuroscience research. The typical duration for activity recording is 5 to 10 minutes per animal. Animals are tested individually to avoid interference from others. The device must be placed in a quiet and undisturbed environment during testing. Both

horizontal and vertical activity can be measured in advanced models. Due to its simplicity, reproducibility, and sensitivity, the actophotometer is widely used in neuropharmacology, toxicology, and behavioural research for rapid screening of CNS-active drugs, dose-response studies, and comparative activity assessments.

The actophotometer works on the principle that when an animal moves inside a closed chamber, it interrupts a beam of light (photoelectric cell), and each interruption is recorded as a count. The instrument contains a series of light beams and corresponding photoelectric sensors arranged in the activity cage. When the animal (usually a rat or mouse) crosses the light beam, the beam is broken, and the photoelectric cell converts this interruption into an electrical signal, which is counted and displayed digitally or mechanically. The total number of beam interruptions in a fixed period reflects the locomotor activity of the animal.

Basic Components of Actophotometer: An actophotometer is an electronic device used to measure the spontaneous motor activity of small animals, typically in pharmacological experiments. The movement of the animal interrupts light beams inside the device, which is then recorded as activity counts. This method is widely used to assess the effect of various drugs on

the central nervous system (CNS), especially stimulants and depressants

Essential components that make up an actophotometer and their respective functions:

1. Enclosure or Activity Cage Structure: The primary part of the actophotometer is the enclosed arena or cage, usually square or circular in shape, designed to restrict the animal's movement to a specific area. Material: Typically made of transparent acrylic plastic or metal, it prevents the animal from escaping while allowing clear observation. Function: It serves as the space where the animal is placed during the experiment. The movements of the animal inside this chamber are the basis for data recording.

2. Infrared (IR) Light Source and Sensors IR Light Beams: Multiple infrared beams are placed across the floor of the enclosure. These beams are invisible to the animal and do not interfere with its natural movement. Photoelectric Sensors: Each beam is aligned with a photo sensor on the opposite side of the enclosure. When the animal moves and breaks a beam, the interruption is detected by the sensor. Function: The system detects the number of beam breaks, which reflects the locomotor activity of the animal. The more active the animal, the more beams are broken.

3. Digital Display or Counter Display Unit: A digital or analogue display unit is connected to the beam-sensor system. It shows the number of beam interruptions over a specified period. Types: Modern actophotometer come with digital LCD displays or computer interfaces for accurate data logging and analysis. Function: This component provides a real-time readout of the animal's activity. The readings help researchers compare activity levels before and after drug administration.

4. Timer Built-in Clock or Timer: A timing mechanism is often integrated with the display unit. Function: It allows the experimenter to set the duration of the activity measurement, usually ranging from a few minutes to an hour. Importance: Consistent time intervals are crucial for comparing control and test results.

5. Electrical Power Supply Power Source: Most actophotometer operate using AC mains power or DC battery systems, depending on the model. Function: Powers the infrared lights, sensors, and display systems. Safety Features: Many systems have fuses or surge protectors to ensure safe operation during experiments.

6. Control Panel Buttons/Switches: The control panel allows the user to turn on/off the device, reset counters, and set timing functions. Advanced Models: Some include data storage, USB ports, or software interfaces for direct computer connectivity. Function: Enables customization of experimental settings and simplifies operation.

7. Calibration and Sensitivity Adjustment Knob Use: Some actophotometer offer the option to calibrate the sensitivity of the IR sensors. Function: Adjusts the threshold at which beam interruptions are recorded, helping to fine tune the system based on the size and type of animal.

Standard Operating Procedure:

1. Ensure the Actophotometer is clean and calibrated.
2. Switch on the instrument and allow it to warm up for 5 minutes.
3. Set the timer (commonly 5 or 10 minutes) for recording activity.
4. Administer the test or control drug to the animal by the recommended route.
5. Allow 30–60 minutes (as per protocol) for drug absorption.



6. Place one animal in the activity chamber and start the timer.
7. Record the activity score (digital count) displayed at the end of the session.
8. Repeat the process for all animals (one at a time).
9. After each test, clean the chamber if soiled and reset the counter.
10. Record and analyse the data to assess CNS depressant or stimulant activity.

Precautions:

Handle animals gently and calmly to avoid stress, ensure infrared sensors are unobstructed during the test, Keep the testing area quiet and distraction-free, use freshly prepared drug solutions, do not place more than one animal at a time.

Animals used:

In experimental pharmacology, the actophotometer is a valuable instrument used to measure the spontaneous locomotor activity of animals. To study such effects, small laboratory animals are used. The selection of the animal depends on factors such as size, handling convenience, availability, cost, and relevance to human physiology.

1. Mice (*Mus musculus*):

Most commonly used species in actophotometer studies. Advantages: Small size, easy to handle, require minimal space and food, highly sensitive to CNS-active drugs, Inbred strains provide genetic consistency.

2. Rats (*Rattus norvegicus*):

Used when larger sample size or more tissue is needed for follow-up biochemical studies. Advantages: Larger body size allows easier administration of drugs and sampling; Behavioral

responses are more stable compared to mice. Some complex tasks or longer-duration activity measurements are more reliable in rats.

Applications of Actophotometer:

1. Measures the general locomotor activity of small laboratory animals.
2. Evaluates the stimulant or depressant action of CNS drugs.
3. Screens newly developed neuroactive or psychotropic agents.
4. Assesses behavioral effects of pharmaceutical compounds.
5. Detects alterations in movement due to neurological disorders.
6. Determines the sedative potential of herbal or synthetic drugs.
7. Assists in analysing the anxiolytic or antidepressant activity.
8. Evaluates neurotoxicity through changes in locomotor behaviour.
9. Helps in studying drug tolerance and physical dependence.
10. Monitors activity levels in genetic or disease model animals.
11. Assesses motor function recovery in neurodegeneration studies.
12. Supports dose-optimization in preclinical pharmacology.



Fig.2 Diagrammatic representation of Actophotometer

CONCLUSION:

The Electroconvulsimeter and actophotometer remain cornerstone instruments in pharmacological and behavioural research due to their reliability, sensitivity, and translational relevance. The Electroconvulsimeter, through the MES model, provides a standardized approach for identifying and validating anticonvulsant agents, while the actophotometer enables precise measurement of drug-induced alterations in locomotor activity. Their combined application not only enhances the understanding of drug mechanisms but also accelerates the development of novel CNS-active therapeutics. By offering reproducible and cost-effective methods, these instruments continue to play a crucial role in bridging experimental findings with clinical potential in the field of neuropharmacology.

REFERENCES

1. Borowicz-Reutt K, Banach M. Chronic Treatment with Oxcarbazepine Attenuates Its Anticonvulsant Effect in the Maximal Electroshock Model in Mice. *International Journal of Molecular Sciences*. 2024;25(12):6751.
2. Pinto LF, Silva LS, João RB, Boldrini V, Cendes F, Yasuda CL. Practices in the prescription of antiseizure medications: Is it time to change? *Arquivos de Neuro-Psiquiatria*. 2024;82:s00431777806.
3. Beydoun A, DuPont S, Zhou D, Matta M, Nagire V, Lagae L. Current role of carbamazepine and oxcarbazepine in the management of epilepsy. *Seizure*. 2020;83:251–63.
4. Billakota S, Devinsky O, Kim KW. Why we urgently need improved epilepsy therapies for adult patients. *Neuropharmacology*. 2020;170:107855.
5. Luszczki JJ, Panasiuk A, Zagaja M, et al. Polygonogram and isobolographic analysis of interactions between various novel antiepileptic drugs in the 6-Hz corneal stimulation-induced seizure model in mice. *PLoS ONE*. 2020;15:e0234070.
6. Borowicz-Reutt K, Banach M. Chronic Treatment with Oxcarbazepine Attenuates Its Anticonvulsant Effect in the Maximal Electroshock Model in Mice. *International Journal of Molecular Sciences*. 2024;25(12):6751.
7. Moezi L, Pirsalami F, Dastgheib M, Oftadehgan S, Purkhosrow A, Sattarinezhad E. Acute and Sub-chronic Anticonvulsant Effects of Edaravone on Seizure Induced by Pentylene-tetrazole or Electroshock in Mice, Nitric Oxide Involvement. *Iran Journal of Medical Sciences*. 2023;48(3):329-340.
8. Zhao J, Liang D, Xie T, Qiang J, Sun Q, Yang L, Wang W. Nicorandil Exerts Anticonvulsant Effects in Pentylenetetrazol-Induced Seizures and Maximal-Electroshock-Induced Seizures by Downregulating Excitability in Hippocampal Pyramidal Neurons. *Neurochemical Research*. 2023;48(9):2701-2713
9. Kaminski S, et al. Enhancement of Glutamate Uptake as Novel Antiseizure Approach: Preclinical Proof of Concept. *Annals of Neurology*. 2025.
10. Kruse SW, Dayton KG, Purnell BS, Rosner JI, Buchanan GF. Effect of monoamine reuptake inhibition and $\alpha 1$ blockade on respiratory arrest and death following electroshock-induced seizures in mice. *Epilepsia*. 2019;60(3):495–507.
11. Klein B, Budiman SN, Homberg J, Verma D, Keijer J, van Schothorst EM. Measuring locomotor activity and behavioral aspects of rodents living in the home-cage. *Frontiers in Behavioral Neuroscience*. 2022;16:877323.



12. Terry S, Gommet C, Kerangueven AC, et al. Activity in group-housed home cages of mice as a novel preclinical biomarker in oncology studies. *Cancers*. 2023;15:194798.
13. Zentrich E, Talbot SR, Bleich A, Häger C. Automated home-cage monitoring during acute experimental colitis in mice. *Frontiers in Neuroscience-Switzerland*. 2021;15:15.
14. Zhang XY, Ye F, Yin ZH, et al. Research status and trends of physical activity on depression or anxiety: a bibliometric analysis. *Frontiers in Neuroscience-Switzerland*. 2024;18:1337739
15. Camilleri MPJ, Bains RS, Williams CKI. Of Mice and Mates: Automated Classification and Modelling of Mouse Behaviour in Groups using a Single Model across Cages.2023.
16. Ren J, Xiao Z, Zhang Y, Yang Y, He L, Yoon E, Bello ST, Chen X, Wu D, Tortorella M, He J. EPIDetect: Video-based convulsive seizure detection in chronic epilepsy mouse model for anti-epilepsy drug screening.2024.

HOW TO CITE: M. Nithya, G. M. Sivakumar, G. Kiruthika, A. Kokila, A. Mohanasundaram, S. Prasanna, N. Rithikadevi, V. Selvaraja, S. R. Shibi, T. Thirisha, Instrumental Approaches in Neuropharmacological Screening: A Review of Insights from Electroconvulsimeter and Actophotometer Studies, *Int. J. of Pharm. Sci.*, 2025, Vol 3, Issue 8, 2670-2677. <https://doi.org/10.5281/zenodo.16940508>

